
Controls – Ranger 3000







| Button | Primary Function (Short Press ¹) | Secondary Function (Long Press ²) | Menu Function (Short Press) |
|---|---|---|--|
| <div style="border: 1px solid black; padding: 2px; text-align: center;">On/Zero Off</div> <p style="text-align: center;">Yes</p> | <ul style="list-style-type: none"> • Turns scale on • Sets display to zero | <ul style="list-style-type: none"> • Turns the scale off | <ul style="list-style-type: none"> • Selects/accepts displayed setting |
| <div style="border: 1px solid black; padding: 2px; text-align: center;">Print Units</div> <p style="text-align: center;">No</p> | <ul style="list-style-type: none"> • Sends the current value to the selected COM ports if AUTOPRINT is set to Off. | <ul style="list-style-type: none"> • Toggles through active weighing units | <ul style="list-style-type: none"> • Advances to the next menu or menu item. • Rejects the current setting on the display and advances to the next available setting |
| <div style="border: 1px solid black; padding: 2px; text-align: center;">Function Mode</div> <p style="text-align: center;">Back</p> | <ul style="list-style-type: none"> • Initiates an application mode | <ul style="list-style-type: none"> • Allows changing the application mode | <ul style="list-style-type: none"> • Returns to previous menu item |
| <div style="border: 1px solid black; padding: 2px; text-align: center;">M+ Menu</div> <p style="text-align: center;">Exit</p> | <ul style="list-style-type: none"> • Accumulates the weight and displays the accumulated information at 0 load. | <ul style="list-style-type: none"> • Enters User menu | <ul style="list-style-type: none"> • Quickly exit User Menu. • Aborts calibration in progress |
| <div style="border: 1px solid black; padding: 2px; text-align: center;">Tare</div> | <ul style="list-style-type: none"> • Enter/clear a tare value. Clears the accumulation when the accumulation information is displayed. | | |

Notes: ¹ Short Press: Press less than 1 second.


² Long Press: Press and hold for more than 2 seconds.

Controls – Ranger Count 3000

| Button | Primary Function (Short Press) | Secondary Function (Long Press) | Menu Function (Short Press) | Library Function (Short Press) |
|---|--|---|--|---|
|  Yes | <ul style="list-style-type: none"> • Turn scale on • Zero the display | <ul style="list-style-type: none"> • Turn scale off. | <ul style="list-style-type: none"> • Accept the current menu or setting. | <ul style="list-style-type: none"> • Accept the current setting. |
|  No | <ul style="list-style-type: none"> • Send the displayed value to the COM port | <ul style="list-style-type: none"> • Change the weighing unit. | <ul style="list-style-type: none"> • Advance to the next menu or setting. • Increment the displayed value. | <ul style="list-style-type: none"> • Advance to the next library or setting. • Increment the displayed value. |
|  Back | <ul style="list-style-type: none"> • Initiate the function of the current application mode. | <ul style="list-style-type: none"> • Switch between Check Weigh, Check Count and Off. | <ul style="list-style-type: none"> • Previous menu • Decrement the displayed value. | <ul style="list-style-type: none"> • Previous library • Decrement the displayed value. |
|  Exit | <ul style="list-style-type: none"> • Accumulate the weight or pieces. | <ul style="list-style-type: none"> • Enter user menu. • View the Audit Trail event counters | <ul style="list-style-type: none"> • Exit the menu. • Abort the calibration | <ul style="list-style-type: none"> • Exit the library. |
|  | <ul style="list-style-type: none"> • ID number input | | | |
|  | <ul style="list-style-type: none"> • Display/Store APW | | | |
|  | <ul style="list-style-type: none"> • Enter numeric values on the display. | | | |
|  | <ul style="list-style-type: none"> • Enter dot (.) on the display. | <ul style="list-style-type: none"> • Switch between scale 1 and scale 2 | | |
|  | <ul style="list-style-type: none"> • Clear last character from the display. | <ul style="list-style-type: none"> • Switch between positive and negative value | | |
|  | <ul style="list-style-type: none"> • Perform a tare operation. | | | |
|  | <ul style="list-style-type: none"> • Display/Start APW | | | |

Controls – Valor 7000

| Button | Primary Function (Short Press ¹) | Secondary Function (Long Press ²) | Menu Function (Short Press) |
|----------|---|---|--|
| Yes | <ul style="list-style-type: none"> • Turns scale on • Sets display to zero | <ul style="list-style-type: none"> • Turns the scale off | <ul style="list-style-type: none"> • Selects/accepts displayed setting |
| No | <ul style="list-style-type: none"> • Sends the current value to the selected COM ports if AUTOPRINT is set to Off. | <ul style="list-style-type: none"> • Toggles through active weighing units | <ul style="list-style-type: none"> • Toggles through available settings |
| Back | <ul style="list-style-type: none"> • Selects function setting | <ul style="list-style-type: none"> • Selects active Mode | <ul style="list-style-type: none"> • Returns to previous settings |
| Exit | <ul style="list-style-type: none"> • Accumulates the weight and displays the accumulated information at 0 load. | <ul style="list-style-type: none"> • Enters User Menu | <ul style="list-style-type: none"> • Quickly exit User Menu |
| | <ul style="list-style-type: none"> • Enter/clear a tare value. Clears the accumulation when the accumulation information is displayed. | | |
| | <ul style="list-style-type: none"> • The IR Sensor can be programmed to act as "touchless" button. | | |

Notes: ¹ Short Press: Press less than 1 second.

² Long Press: Press and hold for more than 2 seconds.

Ohaus Corporation, 7 Campus Drive, Suite 310, Parsippany, New Jersey, 07054, USA www.ohaus.com

Declaration of conformity We, Ohaus Corporation, declare under our sole responsibility that the balance models listed below marked with "CE" – are in conformity with the directives and standards mentioned.

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
For non-automatic weighing instruments used in an Article 1, 2.(a) application, additional metrological marking according to Annex IV of Council directive 2009/23/EC must be attached to the instrument.

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Issued on: 2012-12-04

Ted Xia
President
Ohaus Corporation
Parsippany, NJ USA

Robert Hansen
Compliance Manager
Ohaus Corporation
Parsippany, NJ USA